

## LIFT CHAIR POSITION GUIDE

	Lifted	Seated	Reading	Napping	Relaxing	Sleeping	Trendelenburg	Zero-Gravity
2-Position	<b>*</b>	<b>*</b>	<b>*</b>	×	×	×	×	×
3-Position	<b>*</b>	<b>*</b>	<b>*</b>	<b>*</b>	×	×	×	×
Infinite-Position	<b>*</b>	*	~	<b>*</b>	<b>*</b>	<b>*</b>	×	×
Zero-Gravity	<b>*</b>	<b>*</b>	~	<b>~</b>	~	<b>*</b>	<b>*</b>	~



2-POSITION TV RECLINE The 2-position TV recline is achievable on all Pride® Lift Chairs. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.



**3-POSITION FULL RECLINE** 

The 3-position full recline is available on most Pride® Lift Chair models and puts you in a position ideal for napping or relaxing.



**INFINITE-POSITION** SLEEP RECLINE

The infinite-position sleep recline is perfect for taking a nap – whether quick or extended – or even for a full night's sleep. Infinite-position Pride® Lift Chairs enable the backrest and footrest to move independently from one another so that a variety of positions can be reached.



**TRENDELENBURG** POSITION

The Trendelenburg position, available on the LC-525, SR-525 and LC-576, elevates one's feet higher than his or her head for increased comfort and relaxation.



**ZERO-GRAVITY POSITION** 

The zero-gravity position, which can be reached on infinite-position Pride® Lift Chair models, positions one's thighs and lower legs higher than his or her torso.